
# FORA<sup>®</sup> ADVANCED pro GD40g

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










(t)

## Safety Information

Read the following *Safety Information* thoroughly before using the device.

- Use this device **ONLY** for the intended use described in this manual.
- Do **NOT** use accessories which are not specified by the manufacturer.
- Do **NOT** use the device if it is not working properly or damaged.
- Do **NOT** use the equipment where aerosol sprays are being used or where oxygen is being administered.
- This device does **NOT** serve as a cure for any symptoms or diseases. The data measured is for reference only.
- Before using this device to test blood glucose, read all instructions thoroughly and practice the test. Carry out all the quality control checks as directed.
- Keep the device and testing supplies away from young children. Small items such as the battery cover, batteries, test strips, lancets and vial caps are choking hazards.
- Use of this instrument in a dry environment, especially if synthetic materials are present (synthetic clothing, carpets etc.) may cause damaging static discharges that may cause erroneous results.
- Do NOT use this instrument in close proximity to sources of strong electromagnetic radiation, as these may interfere with the correct operation.
- Proper maintenance and timely calibration with the control solution are essential to the longevity of your device. If you are concerned about the accuracy of the measurement, please contact the place of purchase or customer service representative for assistance.

**KEEP THESE INSTRUCTIONS IN A SAFE PLACE**

## Important Information

- Severe dehydration and excessive water loss may cause readings which are lower than actual values. If you believe you are suffering from severe dehydration, consult a healthcare professional immediately.
- If your blood glucose results are lower or higher than usual, and you do not have symptoms of illness, first repeat the test. If you have symptoms or continue to get results higher or lower than usual, follow the treatment advice of your healthcare professional.
- Use only fresh whole blood sample to test your blood glucose. Using other substances will lead to incorrect results.
- If you are experiencing symptoms that are inconsistent with your blood glucose test results and you have followed all instructions described in this owner's manual, call your healthcare professional.
- We do not recommend using this product on severely hypotensive individuals or patients in shock. Readings which are lower than actual values may occur for individuals experiencing a hyperglycaemic-hyperosmolar state, with or without ketosis. Please consult the healthcare professional before use.
- The measurement unit used for indicating the concentration of blood or plasma glucose can have mg/dL or mmol/L. The approximate calculation rule for conversion of mg/dL in mmol/L is:

mg/dL	Divided by 18	= mmol/L
mmol/L	Times 18	= mg/dL

For example:

1)  $120 \text{ mg/dL} \div 18 = 6.6 \text{ mmol/L}$

2)  $7.2 \text{ mmol/L} \times 18 = 129 \text{ mg/dL}$  approximately.

# Introduction

## Intended Use

This system is intended for use outside the body (in vitro diagnostic use) by people with diabetes at home as an aid to monitoring the effectiveness of diabetes control. It is intended to be used for the quantitative measurement of glucose (sugar) in fresh capillary. It should not be used for the diagnosis of diabetes, or screening for diabetes mellitus.

The test for hematocrit (HCT) as part of the system, is intended for use in the in vitro quantification of packed red blood cell volume fraction in capillary whole blood as an aid in monitoring the status of total volume of red blood cells. The reading of HCT is only for determination of whether the blood test sample is within the acceptable range for the blood glucose system. It should not be used for the diagnosis of anemia or erythrocytosis.

## Test Principle

Your system measures the amount of sugar (glucose) in whole blood. The glucose testing is based on the measurement of electrical current generated by the reaction of glucose with the reagent of the strip. The meter measures the current, calculates the blood glucose level, and displays the result. The strength of the current produced by the reaction depends on the amount of glucose in the blood sample

## Product Overview (a)

- |                       |                        |
|-----------------------|------------------------|
| 1. Test Slot          | 5. Data Port           |
| 2. MAIN Button        | 6. SET Button          |
| 3. Display Screen     | 7. Battery Compartment |
| 4. Test Strip Ejector |                        |

## Screen Display (b)

- |                       |                         |
|-----------------------|-------------------------|
| 1. Blood Drop Symbol  | 9. Measurement Mode     |
| 2. Code               | 10. Measurement Unit    |
| 3. Ketone Warning     | 11. HCT Symbol          |
| 4. Test Result        | 12. HCT Result          |
| 5. Warning Symbol     | 13. Hb Symbol           |
| 6. Low Battery Symbol | 14. Date/Time/Hb Result |
| 7. Day Average        |                         |
| 8. Memory Symbol      |                         |

# Getting Started

## Initial Setup

Before using the device for the first time or after you changed the batteries, do the following initial setup:

### Step 1: Enter the Setting Mode

Open the battery cover and press **SET**. The display turns on.


### Step 2: Configuring the Settings (Date, Time Format, Time, Memory Deletion)

Press **MAIN** repeatedly to adjust the value or enable/disable the setting. Then press **SET** to confirm the setting and switch to another field.

### Note:

- During memory deletion, select **no** to keep all saved results.
- If the device is idle for 3 minutes during the setting mode, it will turn off automatically.

Deletion  
Completed

## Before Testing

### Calibration

You must calibrate the device every time you begin to use a new vial of test strips by setting the meter with the correct code. To ensure test accuracy, make sure the code number displayed on the display screen matches the number printed on the strip vial.

### How to Code Your Meter

1. Insert the code strip into the test slot of the device. Wait for the device to display the code number.

#### Important!

- Make sure the code number on display, code strip, and test strip vial are the same. **(c)**

2. Remove the code strip, the display will show "OFF" indicating the device has finished coding and is ready for blood glucose testing.

### Checking the Code Number

Make sure that the code number displayed on the device matches the number on the test strip vial before you proceed. **(d)**

If it matches, you can proceed with your test. Otherwise, please stop testing and repeat the calibration procedure.

If the problem persists, contact Customer Service for further assistance.

#### Important!

- It is important to make sure that the LCD displayed code is the same as the code on the test strip vial before testing. Failure to do so will get inaccurate results.

## Control Solution Testing

Our Control Solution contains a known amount of glucose that reacts with test strips and is used to ensure your device and test strips are working together correctly.

Do a control solution test when:

- ✓ you first receive the device.
- ✓ at least once a week to routinely check the device and test strips.
- ✓ you begin using a new vial of test strips.
- ✓ you suspect the device or test strips are not working properly.
- ✓ your blood glucose test results are not consistent with how you feel, or if you think the results are not accurate.
- ✓ you have dropped or think you may have damaged the device.

To perform the control solution test, do the following:

1. Insert the test strip into the test slot of the device. Wait for the device to display the test strip and blood drop symbol. **(e)**
2. Press **MAIN** to mark this test as a control solution test.

With "QC" displayed, the device will store your test result in memory. If you press again, the "QC" will disappear and this test is no longer a control solution test. **(f)**

#### Important!

- When doing the control solution test, you have to mark it so that the test result will not be stored in the memory. Otherwise, it will mix up the blood glucose test results with the control solution test results in memory.

## 1. Apply the control solution.

Shake the control solution vial thoroughly before use. Squeeze out a drop and wipe it off, then squeeze another drop and place it on the tip of the vial cap. Hold the device to move the absorbent hole of test strip to touch the drop. Once the confirmation window fills completely, the device will begin counting down. **(g)**

### Note:

- To avoid contaminating the control solution, do not directly apply control solution onto a strip.

## 2. Read and compare the result. **(h)**


After counting down to 0, the test result of control solution will appear on the display. Compare this result with the range printed on the test strip vial and it should fall within this range. If the test result is out of range, read the instructions again and repeat the control solution test.

### Note:

- Do NOT test your blood.
- The control solution range printed on the test strip vial is for control solution use only. It is not a recommended range for your blood glucose level.
- Refer to the **Maintenance** section for important information about your control solutions.
- Out-of-range results:  
If you continue to have test results falling outside the range printed on the test strip vial, the meter and strips may not be working properly.
- There is no HCT and HB display functions while the meter is in QC mode.

# Testing Your Blood Glucose

## Test Strip Appearance

## Inserting a Test Strip

Insert the test strip into its slot.

### Important!

The front side of test strip should face up when inserting test strip. Test results might be wrong if the contact bar is not fully inserted into the test slot.

## Important!

To reduce the chance of infection:

- Never share a lancet or the lancing device.
- Always use a new, sterile lancet. Lancets are for single use only.
- Avoid getting hand lotion, oils, dirt, or debris in or on the lancets and the lancing device.

## Preparing the Lancing Device

1. Remove the cap. **(i)**
2. Insert a new lancet firmly into the white lancet holder cup.
3. Remove the protective disk on the lancet. **(j)**  
Hold the lancet firmly in place and twist off the protective disk.
4. Replace the cap until it snaps or clicks into place.
5. Rotate the dial to set the desired lancing depth. **(k)**
6. Pull the cocking control out until the orange bar appears on the release button window. **(l)**

## Obtaining the Blood Sample

Please follow the suggestions below before obtaining a drop of blood:

- Wash and dry your hands before starting.
- Select the puncture site at fingertips.
- Clean the puncture site using cotton moistened with 70% alcohol and let it air dry.

- Rub the puncture site for about 20 seconds before penetration.

### ❖ Blood from the fingertip

1. Press the lancing device's tip firmly against the lower side of your fingertip.
2. Press the release button to prick your finger, then a click indicates that the puncture is complete. **(m)**

### Note:

- Choose a different spot each time you test. Repeated punctures at the same spot may cause soreness and calluses.
- It is recommended to discard the first drop of blood as it might contain tissue fluid, which may affect the test result.

## Measuring Blood Glucose Test

1. Insert the test strip into its slot.  
Wait for the device to display the test strip and blood drop symbol.
2. Select the appropriate measuring mode by pressing **MAIN**.
  - General Tests (**Gen**) - any time of day without regard to time since last meal.
  - AC (**AC**) - no food intake for at least 8 hours.
  - PC (**PC**) - 2 hours after a meal.
  - QC (**QC**) - testing with the control solution.
3. Obtain a blood sample. (**n**)  
Use the pre-set lancing device to puncture your desired site. After penetration, discard the first drop of blood with a clean tissue or cotton. Gently squeeze the punctured area to obtain another drop of blood. Be careful NOT to smear the blood sample. The volume of blood sample must be at least **1.1** microliter ( $\mu\text{L}$ ) of volume.
4. Apply the blood sample. (**o**)  
Hold the blood drop to touch the absorbent hole of the test strip. Blood will be drawn in and after the confirmation window is completely filled, the device begins counting down. Do not remove your finger until you hear a beep sound.
5. Read your result. (**p**)  
The results of your blood glucose test with hematocrit (HCT) and hemoglobin (Hb) levels will appear after the meter counts down to 0. The results will be stored in meter's memory automatically.

## Disposing Used Test Strip and Lancet

To remove the used test strip, simply press the **Test Strip Ejector** to eject the used test strip. The device will automatically turn off after the test strip is ejected.

To remove the used lancet, remove the lancet from the lancing device after you have finished testing.

Discard your used strip and lancet properly in a puncture resistant container.

### Important!

The used lancet and test strip may be biohazards. Please consult your health care provider for proper disposal complying with your local regulations.

## Reviewing Test Results

Your device stores the 200 most recent test results along with respective dates and times in its memory.

To review all test results, do the following:

1. Press and release **MAIN**. The "**M**" icon appears on the screen.
2. Press **MAIN** to review the test results stored in the device. Press **MAIN** repeatedly to review other test results stored in the device. (**q**) After the last test result, press **MAIN** again and the device will be turned off.

To review the day-average test results, do the following:

1. Press and hold **MAIN** for 3 seconds until the "**DAY AVG**" icon appears. Release **MAIN** and then your 7-day average result measured in general mode will appear on the display.

2. Press **MAIN** to review 14-, 21- and 30- day average results stored in each measuring mode in the order of Gen, AC, then PC. **(r)**

**Note:**

- Press and hold **MAIN** for 5 seconds to exit the memory mode or leave it without any action for 3 minutes. The device will turn off automatically.
- If using the device for the first time, the "----" icon appears when you recall the test results or review the average result. This indicates that there is no test result in the memory.
- Control solution results are **NOT** included in the day average.

## Appendix

### Transferring Data

#### Data transmission via cable

You can use the meter with the RS-232 cable and the Health Care Software System to view test with HCT/Hb levels on your personal computer. To learn more about the Health Care Software System or to obtain the RS-232 cable separately, please contact local customer services or the place of purchase for assistance.

#### Data transmission via cable

##### i. Install the software on your computer

Download Health Care System Software and instruction manual provided on ForaCare Suisse AG, website: <http://www.foracare.ch> to your computer. Follow the instructions to install the software on your computer.

##### ii. Connect the device with your computer

"PC" will appear on the display, indicating that the device is ready to transmit data. **(s)**

##### iii. Transfer data to your computer

Follow the software on-screen instructions to transmit data. Data transmitted will include results with date and time. Remove the cable and the device will automatically turn off.

**Note:**

- When the device is connected to the PC, it is unable to perform any blood glucose test.

## Maintenance

### Changing Battery

When the battery is low, the " " icon appears with E-b. (t)

To change the batteries, do the following:

1. Press the edge of the battery cover and lift it up to remove the cover.
2. Remove the old batteries and replace with two 1.5V AAA size alkaline batteries.
3. Close the battery cover.

When you will hear a "beep" sound, this indicates the batteries are inserted correctly.

#### CAUTION

**RISK OF EXPLOSION IF BATTERY IS REPLACED BY AN INCORRECT TYPE.**

**DISPOSE OF USED BATTERIES ACCORDING TO THE INSTRUCTIONS.**

#### Note:

- Replacing the batteries does not affect the test results stored in memory.
- Keep away these batteries from small children. If swallowed, promptly seek medical assistance.
- Batteries may leak chemicals if unused for a long time. Remove the batteries if you are not going to use the device for an extended period.
- Properly dispose of the used batteries according to your local environmental regulations.

### Caring for Your Device

- To clean the device exterior, wipe it with a cloth moistened with tap water or a mild cleaning agent, then dry the device with a soft dry cloth. Do **NOT** flush with water.
- Do **NOT** use organic solvents to clean the device.

### Device Storage

- Storage condition: -20°C to 60°C (-4°F to 140°F), below 95% relative humidity.
- Always store or transport the device in its original storage case.
- Avoid dropping and heavy impact.
- Avoid direct sunlight and high humidity.

### Caring for Your Test Strips

- Storage condition: 2°C to 32°C (35.6°F to 89.6°F), below 85% relative humidity. Do **NOT** freeze.
- Store your test strips in their original vial only. Do not transfer to other container.
- Store test strip packages in a cool and dry place. Keep away from direct sunlight and heat.
- After removing a test strip from the vial, immediately close the vial cap tightly.
- Touch the test strip with clean and dry hands.
- Use each test strip immediately after removing it from the vial.
- Write the opening date on the strip vial label when you first opened it. Discard remaining test strips after 6 months.
- Do not use test strips beyond the expiry date. This may cause inaccurate results.
- Do not bend, cut, or alter a test strip in any way.

- Keep the strip vial away from children since the cap and the test strip may be a choking hazard. If swallowed, promptly see a doctor for help.

For further information, please refer to the test strip package insert.

### Important Control Solution Information

- Use only our control solutions with your device.
- Do not use the control solution beyond the expiry date or 3 months after first opening. Write the opening date on the control solution vial and discard the remaining solution after 3 months.
- It is recommended that the control solution test be done at room temperature 20°C to 25°C (68°F to 77°F). Make sure your control solution, device, and test strips are at this specified temperature range before testing.
- Shake the vial before use, discard the first drop of control solution, and wipe off the dispenser tip to ensure a pure sample and an accurate result.
- Store the control solution tightly closed at temperatures between 2°C to 30°C (35.6°F to 86°F). Do **NOT** freeze.

### Reference Values

The device provides you with plasma equivalent results.

Time of day	Normal plasma glucose range for people without diabetes (mg/dL)
Fasting and before meal	< 100 mg/dL (5.6 mmol/L)
2 hours after meals	< 140 mg/dL (7.8 mmol/L)

**Source:** American Diabetes Association (2012). *Clinical Practice Recommendations. Diabetes Care, 35 (Supplement 1): S1-100.*

Normal results vary of Hematocrit, but in general are as follows:

- Male: 40.7 to 50.3%
- Female: 36.1 to 44.3%

If the test result of hematocrit is more than 70%, or hemoglobin is more than 23.8 (g/dL), the blood glucose test may be invalid, please redo the test, and consult your doctor if the result is repeatedly invalid.

**Please consult your doctor to determine a target range that works best for you.**

### Symbol Information

Symbol	Referent	Symbol	Referent
	For in vitro diagnostic use		Use within 6 months after first opening
	Do not reuse		Batch code
	Consult instructions for use		Manufacturer
	Keep away from sunlight		Serial number
	Keep dry		Caution, consult accompanying documents
	Temperature limitation		Sterilized using irradiation
	Use by		Do not use if package is damaged

## Troubleshooting

If you follow the recommended action but the problem persists, or error messages other than the ones below appear, please call your local customer service.

### Result Readings

Message	What it Means
	< 10 mg/dL (0.5 mmol/L)
	≥ 240 mg/dL (13.3 mmol/L)
	> 600 mg/dL (33.3mmol/L)

### Error Message

Error Message	Cause	What To Do
	The batteries cannot provide enough power for a test.	Replace the batteries immediately.

	Strip has been used.	Repeat the test with a new strip.
	Problem in operation.	Review the instructions and repeat the test with a new strip. If problem persists, contact the local customer service for assistance.
		
		
		
	You may have removed the strip after applying blood, or insufficient blood volume.	Review the instructions and repeat test with a new test strip.
	HCT>75	
	Wrong code strip is inserted or other coding errors.	Make sure to insert the right code strip for the meter. Check if the code numbers on the code strip and the strip vial label are the same.
	Expired code strip.	Repeat the test with a new code strip.

	Ambient temperature is below the system's operation range.	System operational range is 10°C to 40°C (50°F to 104°F). Repeat the test after the device and test strip have reached the above temperature.
	Ambient temperature is above the system's operation range.	

### Blood Glucose Measurement

Symptom	Cause	What To Do
The device does not display a message after inserting a test strip.	Batteries exhausted.	Replace the batteries.
	Test strip inserted upside down or incompletely.	Insert the test strip with contact bars end first and facing up.
	Defective device or test strips.	Please contact customer service.
The test does not start after applying the sample.	Insufficient blood sample.	Repeat the test using a new test strip with larger volume of blood sample.
	Defective test strip.	Repeat the test with a new test strip.
	Sample applied after the device is automatically turned off.	Repeat the test with a new test strip. Apply sample only when flashing "●" appears on the display.
	Defective device.	Please contact customer service.

The control solution testing result is out of range.	Error in performing the test.	Read instructions thoroughly and repeat the test again.
	Control solution vial was poorly shaken.	Shake the control solution vigorously and repeat the test again.
	Expired or contaminated control solution.	Check the expiration date of the control solution.
	Control solution that is too warm or too cold.	Control solution, device, and test strips should be at room temperature (20°C to 25°C / 68°F to 77°F) before testing.
	Defective test strip.	Repeat the test with a new test strip.
Device malfunction.	Please contact customer service.	

## Specifications

Memory	200 measurement results with respective date and time
Dimensions	110 (L) x 57 (W) x 25 (H) mm
Power Source	Two 1.5V AAA alkaline batteries
Weight	71g
External output	RS-232 Serial Cable
Features	Auto electrode insertion detection Auto sample loading detection Auto reaction time count-down Auto switch-off after 3 minutes without action Temperature warning
Operating Condition	10°C to 40°C (50°F to 104°F), below 85% R.H. (non-condensing)
Storage/Transportation Condition	-20°C to 60°C (-4°F to 140°F), below 95% R.H.
Measurement Units	mg/dL or mmol/L
Measurement Range	10 to 600mg/dL (0.5 to 33.3mmol/L)
Hematocrit range	0~70%
Test Sample	Capillary whole blood
Test Result	Glucose measurements are reported as plasma equivalents

This device has been tested to meet the electrical and safety requirements of: IEC/EN 61010-1, IEC/EN 61010-2-101, EN 61326-1, IEC/EN 61326-2-6, EN 301 489-17, EN 300 328.

## WARRANTY TERMS AND CONDITIONS

With respect to disposable products, ForaCare Suisse warrants to the original purchaser that, at time of delivery, each standard product manufactured by ForaCare Suisse shall be free of defects in material and workmanship and, when used for the purposes and indications described on the labeling, is fit for the purposes and indications described on the labeling. All warranties for a product shall expire as of the product expiration date, or if none, after one (1) year from the date of shipment from ForaCare Suisse. ForaCare Suisse warranty hereunder shall not apply if: (i) a product is not used in accordance with its instructions or if it is used for a purpose not indicated on the labeling; (ii) any repairs, alterations or other work has been performed by Buyer or others on such item, other than work performed with ForaCare Suisse's authorisation and according to its approved procedures; or (iii) the alleged defect is a result of abuse, misuse, improper maintenance, accident or the negligence of any party other than ForaCare Suisse. The warranty set forth herein is conditioned upon proper storage, installation, use and maintenance in accordance with applicable written recommendations of ForaCare Suisse. The warranty furnished hereunder does not extend to damage to items purchased hereunder resulting in whole or in part from the use of components, accessories, parts or supplies not furnished by ForaCare Suisse.

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