FORA<sup>®</sup> ADVANCED pro GD40

312-4272001-023 ver 2.0 2016/02

Blood Glucose Test Strip

## Warnings

- For **in vitro** diagnostic use (for use outside of the body only).
- For single use only.
- Healthcare professionals and other users testing multiple patients with this system should handle everything that has come into contact with human blood carefully to prevent transmitting infectious diseases, including sanitized objects.
- Please read this sheet and your Blood Glucose plus  $\beta$ -Ketone Monitoring System Owner's Manual before you use this test strip. Use only FORA ADVANCED pro GD40 Test Strips with FORA ADVANCED pro GD40 Blood Glucose plus  $\beta$ -Ketone Monitoring System to obtain accurate results, and be covered by the manufacturer's warranty.
- Results may be inaccurate when testing on patients with abnormally low blood pressure or those who are in shock.
- Please do not use FORA ADVANCED pro GD40 Blood Glucose plus  $\beta$ -Ketone Monitoring System on critically ill patients. While the blood glucose result is extremely hyperglycemia (over 600 mg/dL (33.3 mmol/L)), the collection of capillary blood from the approved sample sites is not advised when the peripheral circulation is impaired as the delivery of physiological blood glucose level might not be a true reflection. The following circumstances may apply; severe dehydration as a result of diabetic ketoacidosis or due to stress hyperglycemia, hyperosmolar non-ketotic coma, shock, decompensated heart failure NYHA Class IV or peripheral arterial occlusive disease.
- Keep test strips and lancets away from small children. If swallowed, consult a doctor immediately for advice.

## Intended Use

FORA ADVANCED pro GD40 test strips, when used together with FORA ADVANCED pro GD40 Blood Glucose plus  $\beta$ -Ketone Monitoring System, allow your blood glucose levels to be measured by yourself at home or by healthcare professionals. It uses fresh whole blood samples from the finger. This system is not intended for use in the diagnosis or screening of diabetes mellitus.

Professionals may test with capillary and venous blood sample; home use is limited to capillary whole blood testing. Use only heparin for anticoagulation of whole blood.

## Limitations

- Hematocrit: The hematocrit level is limited to between 0% and 70%. Please ask your healthcare professional if you do not know your hematocrit level.
- Neonatal Use: **This test strip can be used for the testing of newborns.**
- This test strip is used for testing fresh capillary, venous, arterial and neonatal blood.
- Altitude Effects: Altitudes up to 10,742 feet (3,275 m) do not affect test results.
- Please see Appendix: Summary of substances and concentrations in excess of limitation with interference.


## Storage and Handling

**IMPORTANT:** Do not use the test strips if they have expired.

- Test strips expire 6 months after first opening. Write the first opening date on the test strip vial when you first opened it.
- Store the test strips in a cool, dry place between 2°C and 32°C (35.6°F and 89.6°F) and below 85% relative humidity.
- Keep the test strips away from direct sunlight. Do not store the test strips in high humidity.
- Store the test strips in their original vial **ONLY**. Do not transfer them to a new vial or any other containers.
- Do not touch the test strips with wet hands.
- Use each test strip immediately after taking it out of the vial. Close the vial immediately after taking out a strip.
- Keep the vial closed at all times.
- Do not bend, cut, or alter the test strip.

## Testing Your Blood Glucose

**PLEASE WASH AND DRY YOUR HANDS BEFORE PERFORMING ANY TESTS.**

Please refer to your Owner's Manual for more information.

The used lancet and test strip are potentially biohazardous. Please dispose of them carefully according to your local regulations.

## Reading Your Result

Your blood glucose readings deliver plasma equivalent results and are displayed in mg/dL or mmol/L.

## Reference values

Time of day	Normal plasma glucose range for people without diabetes
Fasting and before meal	< 100 mg/dL (5.6 mmol/L)
2 hours after meals	< 140 mg/dL (7.8 mmol/L)

Source: American Diabetes Association (2012). Clinical Practice Recommendations. Diabetes Care, 35 (Supplement 1): S1-100.

**Please consult your doctor to determine a target range that works best for you.**

## Questionable or inconsistent results

- If your test results are unusual or inconsistent with how you are feeling:
  - Make sure the confirmation window of the test strip is completely filled with blood.
  - Check the expiry date of the test strips.
  - Check the performance of your meter and test strip with the control solutions.
- **Please Note:** Unusually high or low blood glucose levels may be symptoms of a serious medical condition. If most of your results are unusually high or low, please contact your healthcare professional.

## Quality Control Testing

Our control solutions contain a known amount of glucose that can react with test strips. You can check the performance of the meter, test strip and your technique by comparing the control solution results with the range printed on the label of test strip vial. Checking regularly can ensure your test results are accurate. Please refer to the Owner's Manual for complete testing instructions.

**IMPORTANT:** The reference range of the control solutions may vary with each new vial of test strips. Make sure you check the range on the label of your current vial.

## Chemical Components

> Glucose dehydrogenase (E. coli) 8% > Enzyme protector 8%  
> Electron shuttle 55% > Non-reactive ingredients 29%

## Additional Information for Healthcare Professionals

Always wear gloves and follow your facility's biohazard control policy and procedures when performing tests involving patient blood samples. Use fresh whole blood samples only. Professionals may use test strips to test capillary and venous whole blood. Sample Size: 1.1  $\mu$ L Reaction Time: 5 seconds  
System Measurement Range: 10 mg/dL to 600 mg/dL  
Hematocrit Range: 0 % to 70 %

## Accuracy

The table below displays how often FORA ADVANCED pro GD40 achieves this target. The chart is based on a study carried out on 100 patients (each patient was tested six times) to see how well FORA ADVANCED pro GD40 performed compared to YSI-2300 reference method results.

For glucose results < 100 mg/dL (5.5 mmol/L), the percentage (number) of FORA ADVANCED pro GD40 results that were within $\pm$ 15 mg/dL (0.83 mmol/L) of the reference method.	100% (192/192)
For glucose results $\geq$ 100 mg/dL (5.5 mmol/L), the percentage (number) of FORA ADVANCED pro GD40 results were within $\pm$ 15% of the reference method.	98.0% (400/408)

## Precision

The CV (%) is less than 5% both in intermediate precision and repeatability.

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