Warnings

- For In vitro diagnostic use (for use outside of the body only).
For single use only.
Healthcare professionals and other users testing multiple patients with this system should handle everything that has comes into contact with human blood carefully to prevent transmitting infectious diseases...

Intended Use

FORA ADVANCED pro GD40 β-Ketone test strips, when used together with FORA ADVANCED pro GD40 Blood Glucose plus β-Ketone Monitoring System Monitoring System, allow your β-Ketone levels to be measured by yourself at home or by healthcare professionals. It uses fresh whole blood samples from the finger. This system is not intended for use in the diagnosis or screening of diabetes mellitus.


Limitations

- Hematocrit: The hematocrit level is limited to between 10% and 70%. Please ask your healthcare professional if you do not know your hematocrit level.
Neonatal Use: This test strip must not be used for the testing of newborns.

Storage and Handling

- IMPORTANT: Do not use the test strips if they have expired.
Store the test strips in a cool, dry place between 4 °C and 30 °C (39.2 °F and 86.0 °F) and below 85% relative humidity.

Testing Your β-Ketone

Please refer to your Owner's Manual for more information. The used lancet and test strip are potentially biohazardous. Please dispose of them carefully according to your local regulations.

Reading Your Result

The β-Ketone readings deliver plasma equivalent results and are displayed in millimoles of β-Ketone per liter of blood (mmol/L). The β-Ketone test measures Beta-Hydroxybutyrate (β-OHB), the most important of the three β-Ketone bodies in the blood. Normally, levels of β-OHB are expected to be less than 0.6 mmol/L.

Please consult your doctor to determine a target range that works best for you. Questionable or inconsistent results If your test results are unusual or inconsistent with how you are feeling... Please Note: Unusually high or low β-Ketone levels may be symptoms of a serious medical condition.

Quality Control Testing

Our control solutions contain a known amount of glucose that can react with test strips. You can check the performance of the meter, test strip and your technique by comparing the control solution results with the range printed on the label of vial or individual foil pack.

IMPORTANT: The reference range of the control solutions may vary with each new test strip. Make sure you check the range on the label of vial or individual foil pack of your current test strip.

Chemical Components

- β-Hydroxybutyrate Dehydrogenase (Pseudomonas sp.) ≥0.5 U
Electron shuttle 35%
Enzyme protector 10%
Non-reactive ingredients 42%

Additional Information for Healthcare Professionals

Always wear gloves and follow your facility's biohazard control policy and procedures when performing tests involving patient blood samples. Accuracy was assessed with 100 patients by comparing whole blood β-OHB results with plasma results obtained using reference laboratory instruments.

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