FORA 6 Blood Glucose Test Strip

312-0000000-050 ver 2.0 2019/02

Warnings

- For **in vitro** diagnostic use (for use outside of the body only).
- For single use only.
- Healthcare professionals and other users testing multiple patients with this system should handle everything that has come into contact with human blood carefully to prevent transmitting infectious diseases, including sanitized objects.
- Please read this sheet and your FORA 6 Plus / FORA 6 Connect Multi-functional meter and FORA 6 Smart blood glucose meter Owner's Manual before you use this test strip.
- Use only FORA 6 blood glucose test strips with FORA 6 Plus / FORA 6 Connect Multi-Functional meter and FORA 6 Smart blood glucose meter to obtain accurate results, and be covered by the manufacturer's warranty.
- Results may be inaccurate when testing on patients with abnormally low blood pressure or those who are in shock.
- Please do not use FORA 6 Plus / FORA 6 Connect Multi-functional monitoring system and FORA 6 Smart Blood Glucose monitoring system on critically ill patients. While the blood glucose result is extremely hyperglycemia (over 600 mg/dL (33.3 mmol/L)), the collection of capillary blood from the approved sample sites is not advised when the peripheral circulation is impaired as the delivery of physiological blood glucose level might not be a true reflection. The following circumstances may apply: severe dehydration as a result of diabetic ketoacidosis or due to stress hyperglycemic, hyperosmolar non-ketotic coma, shock, decompensated heart failure NYHA Class IV or peripheral arterial occlusive disease.
- Keep test strips and lancets away from small children. If swallowed, consult a doctor immediately for advice.

Intended Use

FORA 6 blood glucose test strip when used together with FORA 6 Plus / FORA 6 Connect Multi-functional meter and FORA 6 Smart blood glucose meter allow your blood glucose levels to be measured by yourself at home or by healthcare professionals. It uses fresh whole blood samples from the finger, and the following areas: the palm, forearm and upper arm. This system is not intended for use in the diagnosis or screening of diabetes mellitus.

Professionals may test with capillary and venous blood sample; home use is limited to capillary whole blood testing. Use ONLY heparin for anticoagulation of whole blood samples. Please do NOT use EDTA for anticoagulation.

Limitations

- Hematocrit: The hematocrit level is limited to between 0% and 70%. Please ask your healthcare professional if you do not know your hematocrit level.
- Neonatal Use: **This test strip can be used for the testing of newborns.**
- Please see Appendix: Summary of substances and concentrations in excess of limitation with interference.
- Altitude Effects: Altitudes up to 3,275 m (10,742 ft) do not affect test results.
- Do not test blood glucose during or soon after a xylose absorption test. Xylose in the blood can produce elevated glucose results.

Storage and Handling

- Use each test strip immediately after taking it out of the vial. Close the vial immediately after taking out a strip.
- **IMPORTANT: Do not use the test strips if they have expired.**
- Keep the vial closed at all times.
- Store the test strips in their original vial **ONLY**. Do not transfer them to a new vial or any other containers.
- Do not touch the test strips with wet hands.
- Do not bend, cut, or alter the test strip.
- Store the test strips in a cool, dry place between 2° C and 30° C (35.6° F and 86° F) and below 85% relative humidity.
- Keep the test strips away from direct sunlight. Do not store the test strips in high humidity.

Testing Your Blood Glucose

PLEASE WASH AND DRY YOUR HANDS BEFORE PERFORMING ANY TESTS.

Please refer to your Owner's Manual for more information.

The used lancet and test strip are potentially biohazardous. Please dispose of them carefully according to your local regulations.

Reading Your Result

Your blood glucose readings deliver plasma equivalent results and are displayed in mg/dL or mmol/L.

Reference values

Time of day	Normal plasma glucose range for people without diabetes
Fasting and before meal	< 100 mg/dL (5.6 mmol/L)
2 hours after meals	< 140 mg/dL (7.8 mmol/L)

Source: American Diabetes Association (2012). Clinical Practice Recommendations. Diabetes Care, 35 (Supplement 1): S1-100.

Please consult your doctor to determine a target range that works best for you.

Questionable or inconsistent results

- If your test results are unusual or inconsistent with how you are feeling:
- Make sure the confirmation window of the test strip is completely filled with blood.
 - Check the expiry date of the test strips.
 - Check the performance of your meter and test strip with the control solutions.
- Please Note:** Unusually high or low blood glucose levels may be symptoms of a serious medical condition. If most of your results are unusually high or low, please contact your healthcare professional.

Quality Control Testing

Our control solutions contain a known amount of glucose that can react with test strips. You can check the performance of the meter, test strip and your technique by comparing the control solution results with the range printed on the label of test strip vial. Checking regularly can ensure your test results are accurate. Please refer to the Owner's Manual for complete testing instructions.

IMPORTANT: The reference range of the control solutions may vary with each new test strip. Make sure you check the range on the label of vial or individual foil pack of your current test strip.

Chemical Components

> Glucose dehydrogenase (*E. coli*) 8% > Enzyme protector 8%
> Electron shuttle 55% > Non-reactive ingredients 29%

Additional Information for Healthcare Professionals

Always wear gloves and follow your facility's biohazard control policy and procedures when performing tests involving patient blood samples. Use fresh whole blood samples only. Professionals may use test strips to test capillary and venous whole blood. Sample Size: 0.5 µL
Reaction Time: 5 seconds
System Measurement Range: 10 mg/dL to 600 mg/dL
Hematocrit Range: 0% to 70%

Accuracy

The table below displays how often FORA 6 achieves this target. The chart is based on a study carried out on 100 patients (each patient was tested six times which resulted in 600 test results) to see how well FORA 6 performed compared to Cobas Integra® 400 plus reference method results. [6]

Table 1 Accuracy results for glucose concentration < 100 mg/dL (5.55 mmol/L)

Within ± 5 mg/dL (Within ± 0.28 mmol/L)	Within ± 10 mg/dL (Within ± 0.55 mmol/L)	Within ± 15 mg/dL* (Within ± 0.83 mmol/L)
33.3% (52/156)	73.7% (115/156)	98.7% (154/156)

Table 2 Accuracy results for glucose concentration ≥ 100 mg/dL (5.55 mmol/L)

Within ± 5%	Within ± 10%	Within ± 15%*
69.1% (307/444)	93.9% (417/444)	99.5% (442/444)

Table 3 Accuracy results for glucose concentrations between 40 mg/dL (2.22 mmol/L) to 456 mg/dL (25.3 mmol/L)

Within ± 15 mg/dL or ± 15% (Within ± 0.83 mmol/L or ± 15%)
99.3% (596/600)

Note: * According to the accuracy criteria of EN ISO 15197:2015, 95% of all differences in glucose values (i.e., Cobas Integra® 400 plus reference values minus glucose values of FORA 6) should be within ±15 mg/dL (0.83 mmol/L) for glucose concentration < 100 mg/dL (5.55 mmol/L), and within ±15% for glucose concentration ≥ 100 mg/dL (5.55 mmol/L). When Test Strips results are compared to the reference values, difference values below 100 mg/dL (5.55 mmol/L) are expressed in mg/dL or mmol/L, while those above 100 mg/dL (5.55 mmol/L) are in percentage.

User performance

160 subjects were tested with blood samples taken from the fingertip and the alternative sites, the palm, the forearm and the upper arm. The tables show how well FORA 6 performed compared to Cobas C311 reference method results.

Table 1 Difference distribution for glucose concentration < 100 mg/dL (5.55 mmol/L)

Tested sites	Difference within ± 5 mg/dL	Difference within ± 10 mg/dL	Difference within ± 15 mg/dL
Fingertip	26/43 (60.5%)	37/43 (86.0%)	43/43 (100%)
Palm	27/42 (64.3%)	41/42 (97.6%)	42/42 (100%)
Forearm	31/42 (73.8%)	39/42 (92.9%)	42/42 (100%)
Upper arm	29/42 (69.0%)	38/42 (90.5%)	41/42 (97.6%)

Table 2 Difference distribution for glucose concentration ≥ 100 mg/dL (5.55 mmol/L)

Tested sites	Difference within ± 5%	Difference within ± 10%	Difference within ± 15%
Fingertip	59/117 (50.4%)	102/117 (87.2%)	115/117 (98.3%)
Palm	49/118 (41.5%)	92/118 (78.0%)	118/118 (100%)
Forearm	43/118 (36.4%)	84/118 (71.2%)	115/118 (97.5%)
Upper arm	49/118 (41.5%)	87/118 (73.7%)	116/118 (98.3%)

Precision

In both intermediate precision and repeatability tests, the standard deviation (SD) is within 5 mg/dL (0.28 mmol/L) for each glucose concentration < 100 mg/dL (5.55 mmol/L) and the coefficient of variation (CV) is less than 5% for each glucose concentration ≥ 100 mg/dL (5.55 mmol/L).

For self-testing [REF ACS045](#)

Use only with **FORA 6 Plus / FORA 6 Connect** Multi-functional meter and **FORA 6 Smart** blood glucose meter.

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Substance	Concentration Tested (mg/dL)	Therapeutic / Physiologic Concentration Range (or Upper Limit) (mg/dL)
Pralidoxime Iodide	> 5,0	~ 10 (i.v. Dose 500 mg)

Referenz / Referenz / Référence / Referenz

- [1] McEnroe, J Robert, et al. National Committee for Clinical Laboratory Standards. Interference testing in clinical chemistry; Approved Guideline - 2nd edition. NCCLS : 2005 - EPT-A2, volume 25, number 27.
- [2] EN ISO 15197: 2015 (E): In vitro diagnostic test systems - Requirements for blood-glucose monitoring systems for self-testing in managing diabetes mellitus, Second edition, section 6.4.4 Interference Testing.
- [3] Hardman JG, Limbird LE, Molinoff PB, et al. Goodman & Gilman's the Pharmacological Basis of Therapeutics. 9th ed. New York, NY: McGraw-Hill; 1996: 1570-1571.
- [4] Zaloga GP, Macgregor D. The Critical Care Drug Handbook. 2nd ed. St Louis, MO: Mosby-Year Book; 1997.
- [5] Kost GJ, Vu HT, Lee JH, et al. Multicenter study of oxygen-intensive handheld glucose point-of-care testing in critical care/hospital/ambulatory patients in the United States and Canada. Crit Care Med. 1998;26:581-590.
- [6] Institut für Diabetes Technologie, System accuracy evaluation of FORA 6 Connect Multi-functional Monitoring System based on ISO 15197:2013 & EN ISO 15197:2015. Project No.: IDT-1739(2A)-FS; 2018

Appendix / Anhang / Annexe / Appendice

Symbol / Symbol / Symbole / Symbol	Referent / Bedeutung / Signification / Bedeutung
	<ul style="list-style-type: none"> • In vitro diagnostic medical device • Für die In-Vitro-Diagnostik • Dispositif médical de diagnostic in vitro • Dispositivo medico per uso diagnostico in vitro
	<ul style="list-style-type: none"> • Do not reuse • Nicht wiederverwenden • À usage unique • Non riutilizzare
	<ul style="list-style-type: none"> • Consult instructions for use • Gebrauchsanweisung beachten • Consulter le mode d'emploi • Leggere le istruzioni per l'uso
	<ul style="list-style-type: none"> • Temperature limitation • Temperaturbegrenzung • Plage de températures limites • Limite di temperatura
	<ul style="list-style-type: none"> • Use by • Zu verwenden bis • Date limite d'utilisation • Data di scadenza
	<ul style="list-style-type: none"> • Keep Dry • Vor Feuchtigkeit schützen • Gardez à sec • Tenere in luogo asciutto
	<ul style="list-style-type: none"> • Keep away from sunlight • Vor Sonnenlicht schützen • Tenir à l'écart de la lumière du soleil • Tenere lontano dai raggi del sole
	<ul style="list-style-type: none"> • CE mark • CE-Zeichen • Marquage CE • Marcatura CE
	<ul style="list-style-type: none"> • Batch code • Chargencode • Code du lot • Codice di lotto
	<ul style="list-style-type: none"> • Humidity limitation • Einschränkung bei Luftfeuchtigkeit • Limite d'humidité • Limite di umidità
	<ul style="list-style-type: none"> • Do not use if package is damaged • Bei beschädigter Verpackung nicht verwenden • Ne pas utiliser si l'emballage est endommagé • Non utilizzare se la confezione è danneggiata
	<ul style="list-style-type: none"> • Manufacturer • Hersteller • Fabricant • Produttore
	<ul style="list-style-type: none"> • Model number • Modell-Nr. • Numéro de modèle • Numero di modello
	<ul style="list-style-type: none"> • Quantity • Menge • Quantité • Quantità

ForaCare Suisse AG
Neugasse 55, 9000,
St. Gallen, Switzerland
www.foracare.ch